

DISTRIBUTION: - Customer Specification for Product 15830**Date Added:** 02-JUN-2022**Product No:** 15830

DESCRIPTION: LAMB LEG CHUMP OFF SHANK OFF BONELESS VP

AMENDMENT: Bag amended.

AMENDED: 27-NOV-2023 13:54 by PHERBERT

WEIGHT RANGE: 15 to 27.2

CRITERIA: Ex: H Grade Lambs Only.

TREATMENT: VP

STATE: Frozen

HALAL: Yes

EU: Yes

TARGET YIELD: 15

CARTON CODE: CTNEDGFR

CARTON CODE 2:

END PANEL: LAMB LEG CHUMP OFF BONELESS
LAMSBOUT ZONDER HEUP ZONDER BEEN (DUTCH)
LAMMKEULE OHNE HUFTE OHNE KNOCHEN (GERMAN)
GIGOT D AGNEAU DESOSSE (FRENCH)
LAMM BENFRI STEK (SWEDISH)

DATE REQUIREMENTS: Date Of Expiry - on End Panel - External
Date Of Pack - on End Panel - External
Date Of Slaughter/Production - on End Panel - External
Date Of Expiry - Internal
Date Of Pack - in Full - Internal
Date Of Production. - Internal

EXPIRY: Expiry Code - 1095DY
Use By Type - Slaughter/Production

LABEL DETAILS: Standard End Panel Label (I:\Translations\Cedar Creek Label Formats\End Panel Labels) - External
Special AFFCo Blank Label (I:\Translations\Ticketpro Label Library) - Internal

PACKAGING DETAILS:

R703131 FRZN 300X390 ME	16 EA
S004730 LABELS MPI SEALS	1 EA
S004845 LABEL SYNTHETIC CARTON 90X130	1 EA
S004874 PRODUCT LABEL AFFCO GENERIC	16 EA
S999957 STRAPPING	.001 EA
T601110 CASE AFFCO 170MM (1pce LID)	1 EA

PREPARED FOR: GERMANY
NETHERLANDS
EUROPEAN UNION

PREPARATION: Derived from fresh carcasses to minimum AC&A MIRINZ specifications or equivalent.
Leg - chump off - shank off is prepared from a leg - chump on (item 4800) by the removal of bones, tendons and lymph nodes. The shank is removed at the stifle joint by a cut through the heel muscles of the silverside. The chump is removed by a straight cut approximately 12mm from the hip joint (acetabulum) and at right angles across the primal.
Flank, lymph nodes, cod fat removed.
Tunnel boned
HAM 5070

NOTES: Remove all bloodstains, loose pieces and bevel any excess fat around the cut surfaces. At no point must the external fat exceed 5mm. Free of sawdust and thoroughly cleaned before packing.

HOW PRODUCT IS PACKED:

CUTS PER BAG: one (1).

BAGS PER CARTON: maximum to catchweight.

Carton Strapping: Individual cartons are strapped twice crossways then twice lengthways.

Labelling: Each pack to be labelled with an Affco loaf label generated through

Cedar Creek piece weigh system. (pack is to placed on scale , enter correct standpac number to activate. Apply label to the back of each pack

Must be printed on an Affco generic loaf label.

CUT CODE: PCEEDGFR

CUT CODE 2:

SPECIAL REQUIREMENTS: Maximise yields during aitch bone and femur bone removal. Score marks to be kept to a minimum.

No Title

